





PHOTOPOL PRO



**dentalfarm**

DENTAL EQUIPMENT SINCE 1972


EXCELLENCE.TODAY



**Real-UV spectrum (< 410 nm) VIOLET**

**Near-UV spectrum (> 410 nm) BLUE**

**UV spectrum (320 - 550 nm) both**



A5407 PHOTOPOL PRO	A5407V PHOTOPOL PRO VACUUM
<i>Height: 300mm - Width: 280mm - Depth: 350mm</i>	<i>Height: 300mm - Width: 280mm - Depth: 350mm</i>
<i>Net and Gross Weight: 8,5Kg - 9,5Kg</i>	<i>Net and Gross Weight: 10,5Kg - 12,5Kg</i>
<i>Voltage: 230V - 50/60Hz (different voltage available on request)</i>	<i>Voltage: 230V - 50/60Hz (different voltage available on request)</i>
<i>Absorption: 250W</i>	<i>Absorption: 300W</i>
<i>Rotating plate size: 100x100mm</i>	<i>Rotating plate size: Ø 95mm</i>
<i>Wavelength range: 320-550 nm (peak at 470 nm)</i>	<i>Wavelength range: 320-550 nm (peak at 470 nm)</i>