

5601004 (005) – DISCS FOR STUMP SEPARATING UNIT



Product group	Accessories for stump separating unit
Main features	⇒ Very thin cutting discs ⇒ Cutting circumference lined with industrial diamond ⇒ Available in 2 versions: 5601004 full circumference and 5601005 circumference splitted into 18 equidistant cutting sectors
Warnings	During use, avoid any contact with metal parts susceptible to damage the surface
Dimensions	Diameter 80 mm x 0.25 thick Hole diameter 16 mm
Packaging	Individual