

A3501 DRILLING KIT



Product group	This instrument is not intended as self-operating unit but it can convert a DENTALFARM surveyor into a DRILLING UNIT
Main features	⇒ Double-joint arm to support any micromotor, turbine or hot spatula ⇒ Most suitable to carry out milling operations both on wax and on metal ⇒ Complete set of accessories allowing to carry out: • restrained vertical movement for drilling and boring • restrained horizontal movement for coulisse • locking the articulated joint ⇒ At same time all procedures usually carried out with the manual Surveyor can be easily executed with this instrument
Number of arms	1 – double-joint arm fitted with bivalent head, featuring on one side the micromotor locking pliers and on the other side the slide-guide self-lubricating bush for the pliers holding rod
Overall dimensions	Overall dimensions of the Surveyors are not modified
Net weight	5,2 Kg corresponding to the total weight of the complete unit