

A4701 BUNSEN BURNER with THERMOCOUPLE for NATURAL GAS
A4702 BUNSEN BURNER with THERMOCOUPLE for LIQUID PROPANE GAS (LPG)



Product group	Gas instruments
Main features	⇒ Highest safety assured by the thermocouple valve which stops gas supply in case the flame accidentally extinguishes ⇒ Wide base support for excellent stability ⇒ Possibility to meter the air-gas mixture ⇒ Connection by means of flexible hose with internal diameter 7 or 13 mm
Fuel	A4701 = Natural gas A4702 = Liquid propane gas (lpg) or buthane gas (camping gas)
Working pressure	A4701 = 22 mbar A4702 = 0.1 – 0.5 BAR
Gas consumption	66 l/h
Overall dimensions	150 x 140 x 100h