

A5610 Model Finishing Unit MF-400



Product Group	Equipment for cutting and finishing gypsum or investment
Essential features	1. Interchangeable carbide cutter (available in different shapes) 2. Wide stainless steel working bench 3. Great stability 4. Movable safety shield over the cutter 5. Powerful and fast engine 6. Fit for the finishing of all models made of gypsum or investment 7. Power supply socket and collection sleeve for dust extractor system
Rotating speed	Fix / 2800 rpm
Rotating direction	Clockwise
Cutter dimensions	12mm – 20mm emerging from the working bench
Power tension	230v AC 50/60 Hz
Absorption	280W
Width	210 mm
Depth	280 mm
Height	310 mm
Net weight	12,0 Kg
Gross weight	14,0 Kg