

A5702 CLAMP FOR TWO FLASKS



Product group

Accessories for mobile dentures

Main features

⇒ **Clamp made from machined aluminum casting**

⇒ **Capacity: 2 flasks at the same time**

Dimensions mm

140 x 75 x 280 (max) – max stroke 30 mm

Net weight g

900

A5703 SINGLE CLAMP



Product group

Accessories for mobile prostheses

Main features

⇒ **Clamp made from machined aluminum casting**

⇒ **Capacity: 1 flask only**

Dimensions mm

140 x 75 x 200 (max) – max stroke = 30 mm

Net weight g

700

A5704 STAINLESS STEEL CLAMP



Product group

Accessories for mobile dentures

Main features

- ⇒ **Stainless steel bracket**
- ⇒ **Capacity: 1 or 2 flasks at the same time**
- ⇒ **Ideal for the hydraulic press**

Dimensions mm

150 x 50 x 200 (max) – max stroke 200

Net weight g

700