

A5705 4-ELEMENT FLASK



Product group	Accessories for mobile dentures
Main features	⇒ Machined brass casting ⇒ Totally disassembled, the components are aligned by calibrated plugs ⇒ To be fixed by means of stirrup or put under the press
Dimensions mm	external 135 x 96 x 67h - internal useful 97 x 89 x 55h
net weight g	1500

A5707 SELF-CLAMPING FLASK



Product group	Accessories for mobile dentures
Main features	⇒ Machined brass casting ⇒ Made up of 2 parts aligned by calibrated plugs and closure through screwed rings ⇒ The clamp is not necessary, it can be put under the press
Dimensions mm	external 150 x 94 x 68h - internal useful 99 x 87 x 58h
net weight g	1900