

A6010 Wax dipping pot HOTWAX



Product group

Electrical appliances for hot waxing-up

Main features

- ⇒ Compact in size
- ⇒ Electronic temperature control (microprocessor)
- ⇒ Wide temperature range from 50°C to 120°C, allowing to choose the best temperature according to every type of wax product, preserving its properties.
- ⇒ 3-digit-display showing both the set temperature and ultraflat keyboard
- ⇒ Conversion key (Celsius or Fahrenheit scale)
- ⇒ Hot pot for the immersion of stumps
- ⇒ Protective cover

Line tension

230v A.C.

Absorption

60W

Length

110 mm

Width

140 mm

Hegith

70 mm

Net weight

0,45 Kg (1,0 kg packed)