

PHT (B-R) – LIGHT-CURING RESIN IN PLATES




Product group

Light-curing products

Main features

- ⇒ **Light-curing resin plates ideal for individual or partial impression holders, denture bases, denture test strips, occlusal blocks, for tracking the joint guide pins and for the external copying of clamps or telescopic crowns**
- ⇒ **3 colors available; natural (PHT), blue (PUTB) and pink (PHTR)**
- ⇒ **Size sufficient to achieve a complete upper model**

Shape



Technical data

Thickness: **2 mm**

Flexural strength: **120N / mm²**

Package size

mm 210 x 140 x 80H

Package weight

1500 g

Packaging

Black box with lid – Contents: 50 plates in one color