

A1503 WHISP-R PRO PROFESSIONAL VACUUM CLEANER



Product group	Suction facility for suction of dust and small debris. The equipment cannot be used to suck up liquids.
Essential features	⇒ One-stage, high speed and high power motor ⇒ Three progressive stages of filtration. First stage with paper filter bag, second stage with long-life microfilter cartridge and HEPA filter ⇒ Electromagnetic filter stirrer, which shakes the filter to let dust drop from the filtering surface, thus ensuring longer working performance. ⇒ Automatic or manual mode ⇒ Possibility to connect two units at the same time ⇒ Universal hose adapter for easy connection
Units which can be served by the suction unit	Dry model trimmer, finishing instrument, stump separating unit, sandblaster, working bench
Filtering degree	1μ
Motor power	1000W
Suction volume	45L/s
Noise level	<59 dB (A)
Line tension	230v A.C. 50 Hz (different tensions on request)
Length	240 mm
Width	370 mm
Height	480 mm
Net weight	13,0 Kg